

***In vitro* Angiogenesis Assay: Endothelial Cell Proliferation, EdU Pulse Chase**

1. Add 9,000 cells in 120µl EGM-2(MV) +growth factors –FBS in a 96 well plate. Add 30ul plasma to each well.

120µl of 1:6 Matrigel suspension (45,000 cells/100µl)

2. Allow to recover at 37°C for 18 hours

3. Add 120ul of 2X (20µM) EdU (Component A), incubate 6 hours

4. Remove media and add 100µl of 3.7% **formaldehyde** in PBS. Incubate for 15 minutes at room temperature.

5. Remove fixative and wash each well twice with 300µl PBS

6. Remove wash solution and add 100 µl of 0.1 Triton X-100 in PBS to each well. Incubate for 15 minutes at room temperature

7. Dilute 10X Click-iT EdU buffer in water.

8. Prepare Click-iT reaction cocktail

In order, add

10.3mL 1X Click-iT EdU reaction buffer

480µl CuSO₄ (Component D)

30µl Alexa Fluor azide (Component B)

1.2 ml 1x Click-iT EdU buffer additive

9. Remove Triton X-100, and wash 2 x 300µl PBS.

10. Add 100µl Click-iT reaction cocktail per well

11. Incubate for 30 minutes at room temperature. Protect from light.

12. Remove the reaction cocktail and wash once with 100µl Click-iT reaction rinse buffer (Component F). Then remove the Click-iT reaction rinse buffer, and add 150µl PBS

13. Image.