## Western Blot Technique to Resolve Von Willebrand Factor Multimers and Degradation Fragments from LVAD Plasma

Part I: Electrophoresis and Immunoblotting

<u>Agarose Gel</u> Gel buffer: 0.1%SDS, 0.375 M Tris-Base 1.0% Agarose = 0.10g agarose per 10mL gel buffer Prepare 15ml/cassette

Sample Prep Dilute plasma 1:40 (10µl+390µl) in 1x LDS sample buffer Heat at 70°C for 10 minutes

<u>Electrophoresis</u> Load 15µl per lane Agarose: 60V for 2 hr 30 min Polyacrylamide: 150V for 1 hr 25 min

<u>Transfer</u> iBlot, setting P2 for 8 min 30 sec Block 1 hr in 5% milk-PBS

Probe for vWF Incubate for 1 hr in 10mL of primary anti-human vWF 1:500 in milk-PBS Wash 3 times for 5 minutes in ~20ml PBS Incubate for 1 hr in 10 mL of secondary goat anti-rabbit 1:3,000 in milk-PBS Wash 3 times for 5 minutes in ~20ml PBS

Image Develop with 2mL HRP substrate Image at high resolution at 10 second intervals, and save images as TIFF files

Probe for Albumin Wash 3 times for 5 minutes in ~20ml PBS Incubate for 20 minutes in 10 mL of HRP-conjugated anti-human albumin 1:10,000 in milk-PBS Wash 3 times for 5 minutes in ~20ml PBS Image for <10 seconds Part II: Agarose/Polayacrylamide Transfer iBlot Gel to PVDF Membrane

- 1. Cut all edges from agarose/ polyacrylamide gel and mark top corner
- 2. Open "bottom" portion of iBlot packet
- 3. Place gel onto the top of the iBlot layers
- 4. Wet a blotting sheet with water and place on top of gel
- 5. Place copper "top" sheet of iBlot packet with gel facing down
- 6. White foam sheet goes on the top of the iBlot machine with metal on the top right side
- 7. Close the machine and run at P2 for 8:30 minutes

## 5% Milk buffer blocking

- 1. First block
  - a. 1/500 dilution in 5% milk
    - i. (1)(x) = (1/500)(10ml)
    - ii. .02ml
    - iii. 20µl antihuman + 10ml milk
- 2. Second block
  - a. 1/3000 dilution in 5% milk
    - i. (1)(x) = (1/3000)(10ml)
    - ii. .003ml
    - iii. 3.3µl anti-rabbit + 10ml milk
- 3. Albumin
  - a. 1/10000 dilution in 5% milk
    - i. (1)(x) = (1/10000)(10ml)
    - ii. .001
    - ii. 4µl anti-human albumin + 40ml milk